






COMMON LAB PPE REFERENCE GUIDE

Lab Coats

Type	Physical Characteristics	Applications
<p data-bbox="94 338 245 369">Traditional</p> 	<p data-bbox="488 338 1091 369">100% cotton to 65%/35% polyester/cotton</p> <p data-bbox="488 390 1081 531">General purpose protection. Protects skin and clothing from dirt, dyes, low-hazard chemicals, low-hazard biohazards, and radioactive material.</p> <p data-bbox="488 552 883 583">Lightweight and breathable.</p> <p data-bbox="488 604 1057 674">Not intended for use when fire hazard is present.</p> <p data-bbox="488 695 854 726">Example Source: Medline</p>	<p data-bbox="1120 338 1531 516">Chemistry labs with low flammable solvent use, and no chemical (pyrophorics) or physical hazards likely to pose a risk of fire.</p> <p data-bbox="1120 537 1458 747">Other lab environments without fire hazards, including those with: biological, radiation, physical and animal hazards.</p>
<p data-bbox="94 829 191 861">Barrier</p> 	<p data-bbox="488 829 854 861">99%/1% polyester/carbon</p> <p data-bbox="488 882 1081 978">Provides a good fluid barrier. Sleeves are cuffed, allowing gloves to be worn over sleeves.</p> <p data-bbox="488 999 1091 1104">Does not permit blood or other potentially infectious materials to pass through due to 3-layer construction.</p> <p data-bbox="488 1125 1081 1194">Polyester material is very flammable, will melt and provide fuel source once ignited.</p> <p data-bbox="488 1215 854 1247">Example Source: Medline</p>	<p data-bbox="1120 829 1515 999">Working with human blood, body fluids, tissues, cells or other potentially infectious material which may contain bloodborne pathogens.</p> <p data-bbox="1120 1020 1526 1199">Other lab environments without fire hazards in which cuffed sleeves and lower particle count are advantageous.</p> <p data-bbox="1120 1220 1498 1325">Should never be worn around an open flame or if there is a potential for fire.</p>
<p data-bbox="94 1371 326 1402">Flame Resistant</p> 	<p data-bbox="488 1371 964 1402">100% Nomex (meta-aramid fiber)</p> <p data-bbox="488 1423 1081 1564">Material is intrinsically flame resistant (not chemically treated). Flame resistance properties will not diminish with laundering.</p> <p data-bbox="488 1585 1065 1690">Nomex is resistant to many solvents and corrosives but is a loose weave material; provides superficial fluid barrier only.</p> <p data-bbox="488 1711 854 1743">Example Source: Bulwark</p>	<p data-bbox="1120 1371 1479 1549">Working with water or air reactive chemicals, large volumes of flammable solvents, or potentially explosive chemicals.</p> <p data-bbox="1120 1570 1498 1675">Environments where there is an open flame or the potential for a fire.</p>




COMMON LAB PPE REFERENCE GUIDE

Eyewear

Type	Physical Characteristics	Applications
<p>Safety Glasses</p> 	<p>All safety glasses must comply with ANSI Z87.1-2010 requirements for minimum protection parameters. Polycarbonate lenses block 99.9% of UV light.</p> <p>Safety glasses will provide impact protection from direct trajectory hazards (foreign objects). They will not provide adequate protection for larger volumes of solvent that may bypass lenses.</p> <p>Example Sources: Uvex 3M Fisher</p>	<p>Most common type of eye protection. Lightweight, comfortable and unlikely to fog.</p> <p>Glasses are worn in low hazard environments and situations.</p> <p>Glasses will not protect from dust or liquid running down the face.</p> <p>Glasses allow a minimal gap between eyewear and face (8 mm or less).</p>
<p>Splash Goggles</p> 	<p>All splash goggles must comply with ANSI Z87.1-2010 requirements for minimum protection parameters.</p> <p>Splash goggles provide a better seal around the face, preventing liquid from contacting eyes. Splash goggles are designed with indirect ventilation.</p> <p>Example Sources: Uvex 3M Pyramex</p>	<p>Goggles are worn when there is the potential for splashes of hazardous material or a concentrated presence of foreign material (dust and debris).</p> <p>Goggles should be tested for a snug fit on the face for optimal protection.</p>
<p>Face Shield</p> 	<p>Face shields must comply with ANSI Z87.1-2010 requirements for minimum protection parameters.</p> <p>Face shields are required when there is need for protection of the entire face and throat.</p> <p>Primary eye protection must always be worn under a face shield.</p> <p>Example Source: Pyramex</p>	<p>Face shields are worn in situations where significant splashing or high hazard procedures.</p> <p>Dispensing cryogenics, using strong corrosives or working with potentially explosive material.</p>


COMMON LAB PPE REFERENCE GUIDE

Disposable Gloves

Type	Physical Characteristics	Applications
<p>Latex</p> 	<p>Provides some chemical resistance. Use manufacturer specific chemical resistance chart to verify suitability.</p> <p>Gloves are sensitive to UV radiation and will degrade. Keep covered to prolong life.</p> <p>As with any thin disposable glove, glove should be removed immediately upon chemical contamination.</p> <p>Note that some individuals may have sensitivities or allergies to latex.</p> <p>Example: Microflex – Evolution One</p>	<p>Most commonly used when working with aqueous material or biological hazards.</p> <p>Handling known or potentially infectious material.</p> <p>Working with animals.</p>
<p>Nitrile</p> 	<p>Wider chemical resistance than latex. Some resistance advantages over chloroprene; use manufacturer-specific chemical resistance chart to verify suitability.</p> <p>As with any thin disposable glove, glove should be removed immediately upon chemical contamination.</p> <p>Example: Kimberly-Clark – Purple Nitrile</p>	<p>Chemical, biological or other environment where incidental solvent exposure is possible.</p> <p>Should not be used when handling large amounts of solvent, pyrophoric or explosive material.</p>
<p>Chloroprene (aka Neoprene)</p> 	<p>Wider chemical resistance than latex. Some resistance advantages over nitrile; use manufacturer-specific chemical resistance chart to verify suitability.</p> <p>As with any thin disposable glove, glove should be removed immediately upon chemical contamination.</p> <p>Example: Microflex – NeoPro</p>	<p>Chemical, biological or other environment where incidental solvent exposure is possible.</p> <p>Should not be used when handling large amounts of solvent, pyrophoric or explosive material.</p>




COMMON LAB PPE REFERENCE GUIDE

Reusable Gloves

Type	Physical Characteristics	Applications
<p>Trionic</p> 	<p>Trionic gloves are made by mixing latex, nitrile and neoprene.</p> <p>Generally has the best qualities of each material, though it is still susceptible to UV degradation. Excellent resistance to most acids and low particulate contamination.</p> <p>Review manufacturer chemical resistance data for specific protection information.</p> <p>Example: MAPA – Trionic E-194</p>	<p>Commonly used in electronics manufacturing and clean room environments.</p> <p>Reasonable cost for a reusable glove.</p>
<p>Nitrile</p> 	<p>Thicker version of disposable nitrile gloves; provides significantly better chemical protection.</p> <p>Review manufacturer chemical resistance data for specific protection information.</p> <p>Example: Ansell – AlphaTec (13 mil)</p>	<p>Chemical, biological or other environment where solvent exposure is likely.</p> <p>Gloves are rated to be immersed in solvent.</p> <p>Gloves should be cleaned and allowed to dry after each use.</p>
<p>Chloroprene (aka Neoprene)</p> 	<p>Thicker version of disposable chloroprene gloves, provides significantly better chemical protection.</p> <p>Review manufacturer chemical resistance data for specific protection information.</p> <p>Example: Ansell – Neoprene (18 mil)</p>	<p>Chemical, biological or other environment where solvent exposure is likely.</p> <p>Gloves are rated to be immersed in solvent.</p> <p>Gloves should be cleaned and allowed to dry after each use.</p>

COMMON LAB PPE REFERENCE GUIDE

Reusable Gloves (continued)

Type	Physical Characteristics	Applications
<p>Viton/Butyl</p> 	<p>Developed in collaboration between UCSF and Ansell Ltd. The layered Viton-over-butyl glove provides excellent protection from chloroform and phenol.</p> <p>Combined thickness is as low as 8 mil, more similar to disposable gloves.</p> <p>Review manufacturer chemical resistance data for specific protection information.</p> <p>Example: Ansell – ChemTek Viton/Butyl</p>	<p>Used most frequently during phenol/chloroform extractions.</p> <p>Good chemical protection; can be used in other situations as appropriately indicated by chemical resistance data.</p>
<p>Laminated Film</p> 	<p>Multiple layers of polyethylene (PE) and ethylene vinyl alcohol (EVOH).</p> <p>Good chemical resistance to a wide array of chemicals. Review manufacturer chemical resistance data for specific protection information.</p> <p>Material is slick, may need to use a disposable overglove for dexterity and grip.</p> <p>Examples: Ansell – Barrier North (Honeywell) – Silver Shield</p>	<p>Working with larger volumes of chemicals, hazardous material spills, or where other glove materials are insufficient.</p> <p>Silver shield provides good resistance to methylene chloride (dichloromethane).</p>
<p>Nomex “Flight gloves”</p> 	<p>Made from Kevlar and Nomex, often with leather palm.</p> <p>Provides excellent flame protection for hands, but does not protect from chemicals.</p> <p>Should be worn under an appropriate chemical glove (nitrile or neoprene).</p> <p>Example: Hatch – BNG Flight Glove</p>	<p>Used when working with pyrophoric or potentially explosive material.</p> <p>Should be used in conjunction with Nomex or other flame resistant lab coat.</p>