



The purpose of this fact sheet is to provide examples of biohazard bags that are compliant with California Department of Public Health (CDPH) medical waste regulations. For more information or to add to this list contact EH&S Biosafety kterry@ucsc.edu.

Compliant Biohazard Bags – Primary Containment

Effective January 1, 2015: Biohazard bags must meet CDPH Medical Waste Management Act, California Health and Safety Code Standards including being **marked and certified** by the manufacturer as passing tests for **tear resistance** by the American Society for Testing Materials as per **ASTM D1922** and **impact resistance ASTM D1709** standards.

Size	Manufacturer	Vendor	Quantity	Catalog Number	Price on CruzBuy (as of 2/14/20)
1 Gallon (8" x 12")	VWR	VWR	400/cs	14230-910	66.14
1 Gallon (8" x 12")	VWR	VWR	400/ca	14220-106	156.46
1.6 Gallon (11" x 14")	Medical Action	Amazon	50/box	50-42	
1.6 Gallon (11"x14")	Medegen	Medline	500/cs	50-40	
1.6 Gallon (11"x14")	Medegen	Medline	500/cs	50-42	
3 Gallon (12" x 24")	VWR	VWR	200/cs	14220-090	77.23
5 Gallon (19" x 23")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828244	32.39
5 Gallon (19" x 23")	Medegen	Medline	400/cs	2634	
5 Gallon (19" x 24")	VWR	VWR	200/cs	14230-912	238.70
5 Gallon (19" x 24")	VWR	VWR	200/cs	14220-110	157.22
10 Gallon (23" x 23")	Stericycle	Stericycle	200/cs	KP2323	
10 Gallon (23" x 23")	Medegen	Medline	500/cs	F116	118.19
10 Gallon (23" x 23")	Medegen	Medline	400/cs	F116BX	67.96
10 Gallon (24" x 23")	ULine	ULine	200/cs	S-12984	
10 Gallon (24" x 24")	Southland Medical	Southland Medical	250/cs	SMC122424	

10 Gallon (24" x 24")	Medegen	Medline	150/cs	37-97	
10 Gallon (24" x 24")	Medegen	Medline	500/cs	47-50	
10 Gallon (24" x 30")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828245	49.73
10 Gallon (24" x 30")	VWR	VWR	200/cs	14220-094	277.63
10 Gallon (28" x 23")	Stericycle	Stericycle	500/cs	1523M	
15 Gallon (24" x 30")	Southland Medical	Southland Medical	250/cs	SMC122433	
15 Gallon (24" x 30")	Medegen	Medline	500/cs	F124	127.42
16 Gallon (24" x 32")	Medegen	Medline	250/cs	45-54	
16 Gallon (25" x 35")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828246	64.19
16 Gallon (25" x 35")	Medegen	Medline	200/cs	2635	
16 Gallon (28" x 31")	Medegen	Medline	250/cs	45-30	
16 Gallon (24" x 36")	VWR	VWR	200/cs	14220-096	282.66
20 Gallon (30" x 36")	VWR	VWR	200/cs	11215-816	233.94
20 Gallon (30" x 36")	ULine	ULine	200/cs	S-20849	
20 Gallon (30" x 36")	Medegen	Medline	250/Cs	F145	
22 Gallon (30.5" x 41")	Medegen	Medline	250/cs	F117	
22 Gallon (31" x 41")	Medegen	Medline	100/cs	47-71	
30 Gallon (31" x 38")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828247	79.15
30 Gallon (33" x 39")	ULine	ULine	100/cs	S-12985R	52.23
30 Gallon (39" x 33")	Stericycle	Stericycle	250/cs	2143	
33 Gallon (31" x 43")	Southland Medical	Southland Medical	200/cs	SMC123143	
33 Gallon (31" x 43")	Southland Medical	Southland Medical	100/cs	NON303143	
35 Gallon (33" x 39")	Medegen	Medline	250/cs	F118	
40 Gallon (37" x 48")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828248	60.69
40 Gallon (37" x 50")	Medegen	Medline	150/cs	F130	
40 Gallon (37" x 50")	Medegen	Medline	150/cs	F110	

40 Gallon (38" x 45")	Medegen	Medline	200/cs	F109	
40 Gallon (38"x 45")	Medegen	Medline	100/cs	F169	
40 Gallon (38" x 48")	VWR	VWR	200/cs	14220-098	225.05
40 Gallon (40" x 45")	Stericycle	Stericycle	100/cs	LIN042	
40 Gallon (40" x 45")	Medegen	Medline	250/cs	F149	
40 Gallon (40" x 46")	Medegen	Medline	100/cs	F135	
40 Gallon (40" x 46")	ULine	ULine	100/cs	S-12986R	
40 Gallon (45" x 40")	Stericycle	Stericycle	200/cs	2346	
40 Gallon (40" x 46")	Southland Medical	Southland Medical	50/cs	SMC304046	
44 Gallon (40" x 46")	Southland Medical	Southland Medical	100/cs	SMC123848	
44 Gallon (40" x 46")	Southland Medical	Southland Medical	60/cs	SMC204046	
55 Gallons (40"x 55")	Southland Medical	Southland Medical	50/cs	SMC304055	
55 Gallons (40"x 55")	Medegen	Medline	150/cs	45-45	
60 Gallons (55"x 60")	Medegen	Medline	50/cs	45-75	
23"x8"x41"	Medegen	Medline	100/cs	47-41	
23"x17"x46"	Medegen	Medline	100/cs	47-43	
23"x8"x41"	Medegen	Medline	100/cs	47-45	